**VET EXPERT**  
BASED ON EVIDENCE

# Vet Expert Rapid Test CaniVec-4

**EN PRINCIPLES**  
The Vet Expert CaniVec-4 is a chromatographic immunoassay for the qualitative detection of *Dirofilaria immitis* antigen and antibody against *Ehrlichia canis*, *Borrelia burgdorferi* and *Anaplasma phagocytophilum*/*Anaplasma platys* in canine serum, plasma or whole blood. Letters T and C on the surface of the plate indicate test line (T) and control line (C). The control line (C) is used for procedural control and should always appear if the test procedure is being performed properly and the test reagents of the control line are working. A colored test line (T) will be visible in the result window if the target antigens and/or antibodies are present in the specimen.  
Sensitivity: *Dirofilaria/Ehrlichia/Borrelia/Anaplasma*: 99.9%/98.2%/93.0%/87.5%  
Specificity: *Dirofilaria/Ehrlichia/Borrelia/Anaplasma*: 99.9%/99.9%/98.0%/99.9%

**MATERIALS PROVIDED (5/10 TESTS/KIT)**  
1) 5/10 Vet Expert CaniVec-4 test plates  
2) 1/2 assay diluent bottle (4 ml)  
3) 5/10 anticoagulant tubes  
4) 5/10 disposable pipettes (20 µl)  
5) Instruction for use

**PRECAUTIONS:**  
1) For preliminary diagnostic use only.  
2) For best results, strict adherence to the instructions is required.  
3) The test plate is sensitive to humidity and heat. Perform the test immediately after opening the pouch.  
4) The test components are single-use only.  
5) Do not touch the membrane in the result window.  
6) Do not use the test kit beyond the stated expiration date marked on the package label.  
7) Do not use the test kit if the pouch is damaged.  
8) Do not mix components from different lot numbers.  
9) All samples should be handled as potentially infectious.  
10) Decontaminate and dispose of all samples and materials in accordance with national and local regulations.  
11) Consider the test results as invalid after the specified read out time.

**STORAGE OF THE TEST KIT**  
The test kit may be kept and stored between 2-30°C. The test is stable within the expiration date. DO NOT FREEZE. Do not expose the test to direct sunlight.

**SPECIMEN COLLECTION AND PREPARATION**  
Serum, plasma or whole blood are appropriate for this test.  
[Serum] Collect the whole blood into the collection tube (NOT containing anticoagulants), leave to settle for 30 minutes for blood coagulation and then centrifuge to get serum.  
[Plasma] Collect the whole blood into the collection tube (containing anticoagulants) and then centrifuge to get plasma.  
[Whole blood] Collect the whole blood into the anticoagulant tube (max. vol. 1.5 ml) if anticoagulated whole blood is not immediately tested, it should be refrigerated at 2-8°C and used within 24 hours.  
Serum samples should be stored at 2-8°C. For longer storage, freeze the samples at -20°C or below. Avoid repeated freezing and thawing. Samples containing preservatives will yield inconsistent test results. They must be clarified prior to assaying. Hemolyzed or contaminated samples may give erroneous results.

**PROCEDURE OF THE TEST (LOOK AT THE TOP)**  
All reagents and samples must reach room temperature (15-30°C) before running an assay  
1) Remove the test from the foil pouch and place it on a flat and dry surface.  
2) Using a disposable pipette (included), add 20 µl (marked on the pipette) of whole blood into each sample well, or add 10 µl of serum/plasma into each sample well using micropipette (not included).  
3) Dispense 3 drops of the assay diluent into each sample well.  
4) The sample will change the result window. If it does not appear after 1 min., add one more drop of the assay diluent to the sample well. Read result after 15 min.

**INTERPRETATION OF THE TEST (LOOK AT THE TOP)**  
1) Negative result – only control line (C) appears in the result window.  
2) Positive result – both test (T) and control (C) lines appear in the result window – indicates the presence of target antigens/antibodies. Even a faint line states a positive result.  
**NOTE:** Test strip of Anaplasma Ab can't differentiate between *A.phagocytophilum* and *A.platys*: a positive result indicates presence of antibodies to *A.phagocytophilum* and/or *A.platys*  
3) Invalid result – if the control line (C) does not appear, the result might be considered invalid. The sample should be retested.

**LIMITATIONS OF THE TEST**  
Although the Vet Expert CaniVec-4 test is very accurate, a low incidence of false results can occur. In case of questionable results, other tests should be performed. Definitive clinical diagnosis should not be based on the result of a single test, but must be an effect of extended clinical and laboratory analysis.

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