VET EXPERT
BASED ON EVIDENCE

Vet Expert Rapid Test CPV/CCV/Giardia Ag

EN	PL	DE	FR	ES	PT	IT	CZ	LV	RO	SI
LT	HU	GR	MK	UA	RU	TR	AL	BG		

EN PRINCIPLES

The Vet Expert CPV/CCV/Giardia Ag test is a chromatographic immunoassay for the qualitative detection of canine Parvovirus, Canine Coronavirus and *Giardia duodenalis* (*G. intestinalis*/*G. lamblia*) antigens in canine feces. Letters T and C on the surface of the plate indicate test line (T) and control line (C). The control line (C) is used for procedural control and should always appear if the test procedure is being performed properly and the test reagents of the control line are working. A colored test line (T) will be visible in the result window if the target antigens are present in the specimen.

Sensitivity: CPV 99.9%/CCV 93.1%/Giardia 99.9%

Specificity: CPV 99.9%/CCV 97.5%/Giardia 99.9%

MATERIALS PROVIDED (2/5/20 TEST/KIT)

- 2/5/20 Vet Expert CPV/CCV/Giardia Ag test plates
- 2/5/20 tubes with diluent for CPV/CCV (P)
- 2/5/20 tubes with diluent for Giardia (G)
- 4/10/40 collection swabs
- 4/10/40 disposable pipettes

PRECAUTIONS:

- For veterinary diagnostics only.
- For best results, strict adherence to the instructions is required.
- The test plate is sensitive to humidity and heat. Perform the test immediately after opening the pouch.
- The test components are single-use only.
- Do not touch the membrane in the result window.
- Do not use the test kit beyond the stated expiration date marked on the package label.
- Do not use the test kit if the pouch is damaged.
- Do not mix components from different lot numbers.
- Handle with the sample as potentially infectious.
- Decontaminate and dispose of all samples and materials in accordance with national and local regulations.
- Consider the test results as invalid after the specified read out time.

STORAGE OF THE TEST KIT

The test kit may be kept and stored between 2-30°C. The test is stable within the expiration date marked on the packaging. DO NOT FREEZE. Do not expose the test to direct sunlight.

SPECIMEN COLLECTION AND PREPARATION

The samples from canine feces should be used for this test. The specimens should be tested immediately after collection if samples are not tested immediately. They should be stored at 2-8°C for 24 hours. For longer storage, freeze at -20°C or below. The right amount of feces to collect is shown in the picture at the top.

PROCEDURE OF THE TEST (LOOK AT THE TOP)

- All reagents must be at room temperature (15-30°C) before use.
- Using 2 swabs, collect the sample from canine feces or directly from the rectum.
- Insert one swab into the specimen tube containing 1ml of assay diluent for CPV/CCV (P-C) and the other swab into the specimen tube containing 1ml of assay diluent for Giardia (G). Stir with the swab for at least 10 seconds so the sample and the diluent can mix well. Discard the swabs while squeezing the swab against the tube wall.
- Remove the test plate from the foil pouch, and place it on a flat and dry surface.
- Using the disposable pipettes provided, take the supernatant from the specimen tube. Use different pipettes for assay diluent (P-C) and assay diluent (G). If the feces particles are too big, wait for 1 minute until the large particles are settled, and then take the supernatant of the sample.
- Add four (4) drops into the sample well.
- Read the result after 10 minutes.

INTERPRETATION OF THE TEST (LOOK AT THE TOP)

- Negative result – only control line (C) appears in the result window.
- Positive result – both test line (T) and control line (C) appear in the result window – indicates the presence of target antigens. Even a fade line states a positive result.
- Invalid result – if the control line (C) does not appear, the result might be considered invalid. The sample should be retested.

LIMITATIONS OF THE TEST

Although the Vet Expert CPV/CCV/Giardia Ag test is very accurate, a low incidence of false results can occur in case of weakly positive samples. Other tests conducted on the same samples should be used to confirm the results based on the results of a single test, but must be an effect of extended clinical and laboratory analysis.

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